

ANEXOS (ANNEXES) al documento de depósito de la línea iPSC SYS FiPS120


ANEXO 1: Análisis de microsatélites

ANNEX 1: Microsatellite analysis


STR	Size range	iPSC	Fibroblast
AMEL	104-112 bps	104/109	104/109
D13S317	166-204 bps	175/175	175/175
D16S539	258-310 bps	290/294	290/294
D18S51	279-364 bps	223/223	222/222
D3S1358	94-148 bps	136/136	135/135
D5S818	108-158 bps	130/130	130/130
D7S820	208-252 bps	232/236	232/236
D8S1179	198-252 bps	226/234	226/234
Penta D	365-448 bps	398/408	398/408
TPOX	255-303 bps	269/281	269/281
vWA	120-184 bps	150/150	150/150

ANEXO 2: Test de diferenciación *in vitro*


ANNEX 2: *in vitro* differentiation test



ANEXO 3: Cariotipo
ANNEX 3: Karyotype



ANEXO 4: Test de micoplasma
ANNEX 4: Mycoplasma test



1. IAC+PC
2. IAC
3. C+
4. C-
7. SYS FiPS120-mR6F-2